

# Guided Surgery.

PF 3.5–PF 4.5  
L 6.5–L 12.5

Instructions for use THM61144




# Content

1. At a glance
  - 3 Intended use
  - 3 Restrictions for use
  - 3 Processing
  - 3 Storage
  
2. Cassette and instruments for guided surgery
  - 4 ELEMENT & CONTACT PF 3.5–PF 4.5 / L 6.5–L 12.5, all collar heights
  
3. Planning and preparation
  - 6 Drilling template preparation
  - 6 Drilling
  - 6 Anatomical structures
  
4. Basic preparation (VECTOdrill™)
  - 7 Double-guide principle
  - 7 Drilling
  - 7 Guide sleeve
  
5. Fine preparation (profile drill)
  - 8 Selecting the profile drill
  
6. Implant insertion
  - 9 Taking up implant and engaging for implant insertion
  - 9 Implant insertion: Visual depth control
  - 9 Rotation index: Implant alignment
  - 9 Adapter for handpiece, one-piece
  
7. Drilling sequence
  - 10 General information
  - 10 Example 1, ELEMENT RC, PF 3.5 / L 11.0
  - 10 Example 2, CONTACT RC, PF 4.0 / L 9.5
  
8. Assembly of the guide sleeve
  - 11 Processing
  - 11 Drilling templates
  - 11 Functional testing
  
9. Important
  - 12 Precautionary measures
  
10. Torque values
  - 13 Torque values
  
11. General notes
  - 14 General notes

# 1. At a glance

These instructions apply to all guided surgery PF 3.5–PF 4.5/L 6.5–L 12.5 specific instruments as listed on page 4/5 in this document and/or shown.

 Only use original Thommen Medical drills with a guiding cylinder Ø 4.8 mm in combination with original guide sleeves Ø 4.8 mm.

For more information on guided surgery of other platforms and lengths, please refer to the VECTOdrill™ pilot drilling solution at [www.ifu-tm.com/THM61141](http://www.ifu-tm.com/THM61141).

For more information on the standard surgical procedure including the use of standard instruments, visit [www.ifu-tm.com/THM61141](http://www.ifu-tm.com/THM61141). For more information on the Thommen 3.0 Implant System, visit [www.ifu-tm.com/THM61146](http://www.ifu-tm.com/THM61146).


## Intended use

Thommen Medical instruments are used in combination with the Thommen Medical Dental Implant System in the upper and lower jaw for implant-borne tooth replacement.

## Restrictions of use

See General restrictions of use (see p. 14).

## Cleaning, disinfection and sterilization

 **Warning:**  
Prepare the cassette and instruments before their first use.

### **Single-use products:**

All products, which are supplied in a non-sterile state, must be sterilized before first use, unless stated otherwise. If prosthetic components have not been reprocessed, no cleaning and disinfection is necessary.

### **Multiple-use products:**

All multiple-use products must be cleaned, disinfected and sterilized before first use. An effective cleaning and disinfection are absolutely necessary requirements for an efficient sterilization for re-use.

Steam sterilization is recommended:

- Fractionated vacuum process with at least 3 vacuum steps (with adequate product drying)
- A steam sterilizer compliant with DIN EN 13060 / DIN EN

285 or ANSI AAMI ST79


- According to EN ISO 17665 validated performance assessment
- Maximum sterilization temperature 138°C (280°F), (plus tolerance in compliance with DIN EN ISO 17665)

Sterilization time, i.e. exposure time at the sterilization temperature, is at least 4 minutes at 132°C (270°F) or 18 minutes at 134°C (273°F) for prion inactivation (not relevant for USA).

For further instructions on the sterilization of prosthetic components, please refer to the valid Thommen Medical processing instructions ([www.ifu-tm.com/THM61131](http://www.ifu-tm.com/THM61131)).

## Storage

The cassette (Art. No. 1.04.020) must be protected from exposure to strong light or high heat during storage.






# 2. Cassette and instruments for guided surgery

## ELEMENT & CONTACT PF 3.5–PF 4.5 / L 6.5–L 12.5, ALL COLLAR HEIGHTS

### Cassette for guided surgery 1.04.020


**Warning:**  
Special slots are provided for loading the cassette with instruments. Pay attention to the dental shaft and graphic template.

#### Example:




#### Basic preparation Initial drill, guided

	L 4.0
Length mm	28.5
Ø 2.0	3.04.110



#### Basic preparation VECTOdrill™ pilot drill and drill, guided

	L 6.5	L 9.5	L 12.5	L 8.0	L 11.0
Length mm	31.5	34.5	37.5	33.0	36.0
Ø 2.0	3.04.111	3.04.113	3.04.115	3.04.112	3.04.114
Ø 2.8	3.04.121	3.04.123	3.04.125	3.04.122	3.04.124
Ø 3.5	3.04.131	3.04.133	3.04.135	3.04.132	3.04.134



#### Fine preparation profile drill, guided

	CONTACT MC	CONTACT RC	ELEMENT MC
Length mm	34.0	32.5	30.5
Ø 3.5	3.04.171	3.04.161	3.04.151
Ø 4.0	3.04.172	3.04.162	3.04.152
Ø 4.5	3.04.173	3.04.163	3.04.153



### MONO torque ratchet\*

Titanium alloy

3.03.160



### MONO guide key\*

Stainless steel

3.03.203



### Guide sleeve

Stainless steel


Length mm	5.0
Ø 4.8	3.04.080 Q4



### Seating instrument for guide sleeve

Stainless steel


Ø 4.8	3.04.100
-------	----------



### Optional instruments\*

In addition to the clearly designated slots, additional slots are available for other instruments.

\* Standard instruments (see [www.ifu-tm.com/THM611141](http://www.ifu-tm.com/THM611141) and [www.ifu-tm.com/THM311111](http://www.ifu-tm.com/THM311111))



### MONO insertion device\*

Stainless steel/PEEK

3.03.162	3.03.163
----------	----------

### MONO screwdriver\*

Stainless steel/PEEK

3.03.165	3.03.166	3.03.167
----------	----------	----------




Max. 20 use cycles (operation, cleaning, sterilization) for cutting instruments. Further use cycles no longer guarantee the function of the product.



Cooling during drilling, approx. 5° C/ 41 °F.


15 rpm



### Implant insertion adapter, guided

Stainless steel

	Adapter, guided	Screwdriver short	Screwdriver extra short
Length mm	25.5	22.0	17.0
	3.04.090	3.03.501	3.03.500



### Adapter for contra angle, one-piece

Stainless steel

	Ø 3.5	Ø 4.0	Ø 4.5-6.0
Length mm			
24.0	3.03.248	3.03.251	3.03.241

# 3. Planning and preparation

## Drilling template preparation

Before the operation, the drilling template must be checked for correct function (correct fit, stability, space) in the patient's mouth. The alignment of the guide sleeves must correspond to the preoperative planning. The respective manufacturer of the drilling template guarantees the compatibility of Thommen Medical guide sleeves with the drilling template. See p. 11 for assembly of the sleeves.

**Planning software:** You can find information on the compatible planning software solutions at [www.thommenmedical.com](http://www.thommenmedical.com)

## Drilling

**Warning:**


All VECTOdrill™ drills are 0.5 mm longer apically than the specified length of the respective Thommen Medical implants. In order to avoid complications, this must be taken into account when choosing and positioning of the implant, particularly in proximity to anatomical structures.

**Note:** If the guided surgery cannot be performed as planned, the guided instruments are also suitable for the standard procedure:

- Use the depth marks for depth control as shown on the graphic insert
- Start with the VECTOdrill™ pilot drill, guided; do not use the initial drill, guided
- Sequential processing with VECTOdrill™ drills and pilot drills, guided, if necessary
- Use the adapter, guided for implant insertion

**Warning:**


Check the drilling depth of the implant bed after the first pilot drill step using the respective laser markings on the drill or using the depth gauge.





## Anatomical structures

The mesiodistal position of the implants can be easily estimated by using a periodontal probe placed vestibularly or determined with a gauge. The required minimum gap width for the sleeve and minimum height for the instruments must be considered.

In case of tight spaces or guided pilot drilling, use the guide sleeve for VECTOdrill™ Pilot Drill Ø 2.0 (Art. No. 3.03.141).



# 4. Basic preparation (VECTOdrill™)




## Double-guide principle

All VECTOdrill™ drills feature a stepped guide tip corresponding to the diameter of the preceding drill. This axis guidance prevents slippage of the drill and ensures a precisely shaped implant bed.



### Warning:

No drilling sequences may be skipped to ensure the double-guide principle.




## Drilling



### Warning:

All holes must be drilled while constantly cooling the exterior with physiological, sterile, cooled saline solution (approx. 5°C/41°F). Recommended rotation speeds must be adhered to to avoid overheating the bone tissue and possible instrument fractures. Lateral forces must be avoided to prevent retention of the guide sleeve in the drilling template. Do not start drilling until the drill is positioned on the bone. Ensure continuous rotation of the drill during extraction.




## Guide sleeve





### Warning:

The guide sleeves are intended only for single use. Multiple use does not guarantee correct functioning of the sleeve.




# 5. Fine preparation (profile drill)





## Selecting the profile drill

The profile drill is selected according to implant type and platform size as specified on the graphic insert.






# 6. Implant insertion



## Take up the implant and engage for implant insertion


**Warning:**

You have to **CLICK!** to ensure correct depth control during implant insertion. **Also check this visually.**




Use the same procedure for the standard implant packaging.

## Implant insertion: Visual depth control




## Rotation index: Implant alignment




## Adapter for handpiece, one-piece

The one-piece adapter for hand piece engages directly into the implant hexagon. It can be used after initially positioning the implant in the bone after the insertion aid has been removed.




# 7. Drilling sequence

## General information




- I Always start with the initial drill
- II Select the length
- III Select the implant type and PF
- IV Identify the last VECTOdill™ drill and follow the lines on the graphic insert
- V Use a profile drill, if necessary
- VI Use the adapter for implant insertion

## Example 1, ELEMENT RC, PF 3.5/ L 11.0



ELEMENT RC: No profile drill required  
(See p. 8).

## Example 2, CONTACT RC, PF 4.0 / L 9.5



# 8. Assembly of the guide sleeve

## Reprocessing

The guide sleeve is suitable for the sterilization procedure described in the processing instructions ([www.ifu-tm.com/THM61131](http://www.ifu-tm.com/THM61131)).



### Warning:

Process the drilling template before use in accordance with the manufacturer's instructions. Make sure that the guide sleeves are securely and correctly fixed in the drilling template before use.

## Drilling templates

The diameter of printed sleeve cavities may vary depending on the 3D printing system used. If diameter corrections are necessary, you can refer to the design settings, if available. You can find recommended printing systems at [www.thommenmedical.com](http://www.thommenmedical.com).



### Warning:

The seating instrument MUST be used for assembly of the sleeve and functional testing. The seating instrument is intended for multiple use. The product must be replaced as soon as it shows signs of wear and tear and/or damage.

Handle for dental technicians

3.03.250


Seating instrument for guide sleeves

Ø 4.8 mm


3.04.100

Guide sleeves


Ø 4.8 mm/L 5.0 mm 3.04.080 Q4



1







2



3

## Functional testing



### Warning:

- The sleeve must be seated securely
- The seating instrument must fit and slide
- The sleeve edge must be correctly positioned vertically

# 9. Important











## Precautionary measures

---



- Process the cassette and instruments before first use.
- You can find processing instructions at [www.ifu-tm.com/THM61131](http://www.ifu-tm.com/THM61131).
- Max. 20 use cycles (operation, cleaning, sterilization) for cutting instruments. Further use cycles no longer guarantee the function of the product.
- Special slots are provided for loading the cassette with instruments. Pay attention to the dental shaft and graphic template.  
Process the drilling template before use in accordance with the manufacturer's instructions.  
Make sure that the guide sleeves are securely and correctly fixed in the drilling template before use.
- Defective or worn instruments must be replaced before the operation.
- Due to the function and construction of the drills, all VECTOdrill drills are 0.5 mm longer apically than the specified length of the respective Thommen Medical implants. In order to avoid complications, this must be taken into account when choosing and positioning of the implant, particularly in proximity to anatomical structures.  
Only use original Thommen Medical drills with a guiding cylinder B 4.8 mm in combination with original guide sleeves Ø 4.8 mm.
- Do not start drilling until the drill is positioned on the bone. Ensure continuous rotation of the drill during extraction.
- Recommended rotation speeds must be adhered to in order to avoid overheating the bone tissue and possible instrument fractures.
- All holes must be drilled while constantly cooling the exterior with physiological, sterile, cooled saline solution (approx. 5°C/41°F). Regularly remove the bone chips to ensure ideal drilling performance.
- Check the drilling depth of the implant bed after the first pilot drill step using the respective laser markings on the drill or using the depth gauge.  
However, when implants are placed in extremely hard bone (such as, for example, in the case of a strongly atrophied edentulous mandible), thread-tapping can be indicated. You can find information on this at [www.ifu-tm.com/THM61141](http://www.ifu-tm.com/THM61141).
- You have to CLICK! to ensure correct depth control during implant insertion. Also check this visually.
- The guide sleeves are intended only for single use. Multiple use does not guarantee correct functioning of the sleeve.
- The seating instrument **MUST** be used for assembly of the sleeve and functional testing. The seating instrument is intended for multiple use. The product must be replaced as soon as it shows signs of wear and tear and/or damage
- Functional testing: The sleeve must be fixed, the seating instrument must fit and slide and the sleeve edge must be correctly positioned vertically.












# 10. Torque values

ABUTMENT/SCREW	TORQUE VALUES IN Ncm				
	10	15	20	25	30
Healing cap Gingiva former 	PF 3.0				
	PF 3.5				
	PF 4.0				
	PF 4.5				
	PF 5.0				
	PF 6.0				
Abutment screws 		PF 3.0			
		PF 3.5			
				PF 4.0	
				PF 4.5	
				PF 5.0	
Dynamic abutment screw 		PF 3.5			
				PF 4.0	
				PF 4.5	
				PF 5.0	
Abutment screw in combination with ART abutment 		PF 3.5			
			PF 4.0		
			PF 4.5		
			PF 5.0		
VARIOmulti abutment Novaloc® abutment 		PF 3.5			
				PF 4.0	
				PF 4.5	
				PF 5.0	
VARIOmulti protective cap 	PF 3.5				
	PF 4.0				
	PF 4.5				
	PF 5.0				
	PF 6.0				
Occlusal screw (VARIOmulti, VARIO 17°, bar) Occlusal closure screw, CAD/CAM bar 		PF 3.5			
		PF 4.0			
		PF 4.5			
		PF 5.0			
Retentive anchor 				PF 4.0	
				PF 4.5	
				PF 5.0	
ZEST LOCATOR® abutment 			PF 3.5		
					PF 4.0
					PF 4.5
					PF 5.0
ZEST LOCATOR® abutment and collar for VARIOmulti 			PF 4.0		
			PF 4.5		
			PF 5.0		
			PF 6.0		

New abutment screws must always be used for permanent insertion of the supraconstruction in the mouth.

# 11. General notes

## THOMMEN IMPLANT SYSTEM

	Manufacturer: Thommen Medical AG Neckarsulmstrasse 28 2540 Grenchen, Switzerland www.thommenmedical.com
	Batch code
	Use by date
	Date of manufacture
	Sterilized using irradiation
	Authorized representative
	Temperature limitation
	Do not re-use
	Non-sterile
	Caution
	Article number
	Conformity symbol as specified by EU Directive MDD 93/42/EEC
	Consult instructions for use
	Do not re-sterilize
	Do not use if package is damaged
	Atmospheric pressure limitation
	Manufacturer
	Keep away from sunlight
	May only be sold to and prescribed by physicians (USA)
	Medical device
	Single product code

### COLORED WARNING STICKER

Application was changed – follow the directions in the corresponding instructions for use.

#### NEW HANDLING

New design – the application has not been changed.

#### NEW DESIGN

**PRODUCT INFORMATION** The information in this document describes the application of the Thommen Medical implant system. This information is available in electronic form online at: [www.ifu-tm.com](http://www.ifu-tm.com). The responsible country representative or distributor for Thommen Medical AG is available to provide technical advice.

**COLOR CODE** Each implant platform diameter has a color code, which can be found on all implant packagings, on the impression items and on most diameter-specific instruments.

### TRACEABILITY

In order to ensure the traceability of the implantable products as well as the manufacturer, product type and product dimensions for a later prosthetic re-restoration, each product package comes with three patient labels. These labels should be used in the practice for documentation and for the implant passport.

Brown	=	PF 3.0
Yellow	=	PF 3.5
Green	=	PF 4.0
Blue	=	PF 4.5
Grey	=	PF 5.0
Purple	=	PF 6.0

**AVAILABILITY** Not all of the Thommen Medical products mentioned in these instructions for use are available in all countries. The responsible country representative or distributor of Thommen Medical AG informs about availability of Thommen Medical products for the country in question.

**GENERAL RESTRICTIONS OF USE** Restorations with cantilevers to individual implants are not recommended. Individual restorations with angled abutments should not be used in regions with high mechanical stress. For implants with a small diameter (PF 3.0 and 3.5), the prosthetic restoration should be constructed in such a way that large bending moment does not occur.

**CONTRAINDICATION** The Thommen Medical products may not be used on patients who are known to have allergies to the corresponding materials.

**POSSIBLE COMPLICATIONS** A stressed loading of the implant or abutment over and above its functional capacity can lead to excessive bone loss or fracture of the implant or restoration. The clinician must supervise the occlusion and functional loading of the prosthetic supraconstruction very carefully.

**SIDE EFFECTS** The patient should be informed about the possible side effects, interactions, precautionary measures and complications associated with Thommen Medical products. Potential complications can occur immediately after insertion of dental implants:

Temporary symptoms: swelling, difficulties with speaking, gum inflammations, pain.

Longer lasting symptoms: chronic pain connected with the dental implant, localized or systemic infections, dysesthesia, loss of alveolar ridge (upper and lower jaw), oroantral or oronasal fistulas, irreversible damage to neighboring teeth, esthetic problems, nerve damage, hyperplasia.

**WARNINGS** All Thommen Medical products that come into effect inside the oral cavity must be protected against aspiration. Thommen Medical products have not been tested for safety and compatibility in an MR environment. Thommen Medical products have not been tested for heating or migration in the MR environment. The safety of Thommen Medical products in the MR environment is unknown. Magnetic resonance tomographic examinations of patients, who have been treated with Thommen Medical products, may result in patient injuries.

**RESPONSIBILITY/LIABILITY** As a part of an overall scheme, Thommen Medical products may be used only with the related original components and instruments in accordance with the instructions for use provided by Thommen Medical. The use of non-system parts may compromise the performance of Thommen Medical products and lead to failures. Users must have appropriate knowledge and information about the handling of Thommen Medical products in order to use the products safely and correctly. The user is obliged to use the Thommen Medical products according to the instructions for use and to check whether the product is suitable for the individual patient situation. The use of Thommen Medical products is the responsibility of the user, as such, beyond the control of Thommen Medical AG. We refuse to accept any responsibility or liability for any damage due to incorrect utilization of the product. Products labeled «Do not re-use» may not be refurbished and/or reused. The refurbishment and/or reuse of these products can affect their function (e.g. fitting and/or cutting properties) as well as their safe use (e.g. risk of infection, disease transmission, fading of the laser or color marks, corrosion). Detailed information

about the possible consequences, which may result from incorrect use, is available from the responsible country representative or distributor of Thommen Medical AG. All serious incidents which have occurred in connection with the product must be reported to the manufacturer and the competent authority of the Member State, in which the user is resident.

**GUARANTEE** The comprehensive guarantees can be found in the country-specific guarantee leaflets.

**TRANSPORT AND STORAGE** Please note the specifications on the labels and instructions for use regarding transportation, storage and handling. If the packaging is damaged, the products must not be used; a visual inspection is necessary. Under no circumstances may Thommen Medical products be used beyond the expiry date, as proper functioning or sterility of sterile packaged products cannot be guaranteed by the manufacturer.

**APPLICATION** The following descriptions are not intended as comprehensive for the immediate use of the Thommen Medical Implant System. Training by a specialist experienced in the use of this system is recommended.

**GUARANTEE OF STERILITY** In general, products of the Thommen Implant System supplied in sterile packaging must not be re-sterilized. Sterile-packed products, whose packaging is damaged, must not be used under any circumstances. Sterile-supplied products, which have not been used for the surgical operation, whose packaging has been opened are considered as having been used and must not be used thereafter. In the event of re-sterilization, proper function and the sterility cannot be guaranteed by the manufacturer. The products intended for single use must never be re-processed, sterilized or reused and must be disposed of safely and properly after use in compliance with all applicable legal and regulatory requirements. Reusable products must be reprocessed according to the instructions for use and, if used on patients, sterilized. They must be checked for their integrity before each use. Any damage (for example, scratches, cracks, nicks, dents), as well as bent parts, means that they must not be used any longer. The number of reprocessing cycles is limited and must be monitored. If the number of cycles is exceeded, proper function and sterility of the product are not guaranteed by the manufacturer anymore.

**DISPOSAL** In the case of cutting products, there is always a risk of injury, therefore the products must be disposed of safely and properly after use, observing all applicable legal and regulatory requirements. For products and their accessories, which have been used on a patient, there is a risk of an infection. Our products are designed and produced so that they can be disposed of safely and correctly after use in compliance with all valid legal and regulatory requirements.

**COPYRIGHT®/REGISTERED®** SPI®, INICELL®, and APLIQUIQ® are registered trademarks of Thommen Medical AG. Publication or reproduction is permitted only with the written consent of Thommen Medical AG. LOCATOR® is a registered trademark of Zest Anchors Inc., CA, USA. Novaloc® is a registered trademark of Valoc AG, Rheinfelden, Switzerland. Dynamic Tibase® is a registered trademark of Talladium España S.L., Lleida, Spain.

**VALIDITY®** Thommen Medical AG. All previous versions lose their validity with the publication of this instruction for use.



## HEADQUARTERS

Thommen Medical AG  
Neckarsulmstrasse 28  
2540 Grenchen | Switzerland  
Tel. +41 61 965 90 20  
Fax +41 61 965 90 21  
info@thommenmedical.com

## SUBSIDIARIES/NATIONAL DISTRIBUTORS

### AUSTRIA

Thommen Medical Austria GmbH  
Mühlgasse 3  
2322 Zwölfaxing | Austria  
Tel. +43 660 2011953  
info@thommenmedical.at

### BENELUX

Thommen Medical Benelux B.V.  
Dierenriem 1  
3738 TP Maartensdijk | Netherlands  
Tel. +31 30 68 68 468  
Info.benelux@thommenmedical.nl

### CHINA

Shanghai Yujing Trading Co., Ltd.  
Room G | Floor 15th | Plaza JiaFa | No.1  
Lane 129 | DaTian Road | JingAn District  
Shanghai | China  
Tel. +86 21 62723077  
Fax +86 21 62175264

### CZECH REPUBLIC

C. Witt Dental spol. s r.o.  
Cihlářská 643/19  
602 00 Brno  
Tel. +420 739 043 449  
helena.novak@cwittdental.cz

### FINLAND

Vector Laboratories Oy  
Engelinaukio 8 B  
00150 Helsinki | Finland  
Tel. +358 400 940 700  
labs@vektor.fi

### FRANCE

Thommen Medical France  
10 avenue Gabriel Pierré  
77680 Roissy-en-Brie | France  
Tel. +33 1 83 64 06 35  
Fax +33 3 89 33 52 53  
infos@thommenmedical.fr

### GERMANY

Thommen Medical Deutschland GmbH  
Am Rathaus 2  
79576 Weil am Rhein | Germany  
Tel. +49 7621 422 58 30  
Fax +49 7621 422 58 41  
info@thommenmedical.de

### HONG KONG

Shengyuan (Hong Kong) Int. Trade Co. Ltd.  
Level 13, 68 Yee Wo Street  
Causeway Bay | Hong Kong  
Tel. +852 530 876 41

### ITALY

Dental Trey S.r.l.  
Via Partisani, 3  
47016 Fiumana | Predappio (FC) | Italy  
Tel. +39 0543 929111  
Fax +39 0543 940659  
implantologia@dentaltrey.it  
www.dentaltrey.it

### JAPAN

J. Morita Corporation  
3-33-18, Tarumi-cho  
Suita | Osaka 564-8650 | Japan  
Tel. +81 6 6384 6921  
Fax +81 6 6384 6746  
www.morita.com

### LITHUANIA/LATVIA

ČERNIKIS MEDICAL PROJECTS, UAB  
Šiaurės prospektas 5B, Kaunas  
Lithuania LT-49191  
Tel. +370 37 201072  
Mobile +370 65 771550  
info@cmp.lt  
www.cmp.lt

### MIDDLE EAST

Star Science International GmbH  
Jupiterstrasse 57  
3015 Bern | Switzerland  
Tel. +41 31 941 07 31  
Fax +41 31 941 07 33  
star.science@bluewin.ch

### NORWAY

Novus Dental AS  
Johannes Bruns gate 5  
0452 Oslo | Norway  
Tel. +47 951 07 007  
post@novusdental.no  
www.novusdental.no

### POLAND

C.WITT DENTAL Sp. z o. o.  
Ul. Granitowa 10  
87-100 Toruń | NIP 951-15-08-371 | Poland  
Tel. +48 56 623 61 23  
biuro@cwittdental.pl  
www.cwittdental.pl

### REPUBLIC OF CROATIA

Futura Dental d.o.o.  
Kralja Zvonimira 108  
10 000 Zagreb | Republic of Croatia  
Tel. +385 91 6814 860  
info@futura-dental.hr  
www.futura-dental.hr

### RUSSIAN FEDERATION

CIS – JSC Geosoft  
Build. 14, Ap. 16, 3-ya Mytishchinskaya ul.  
Moscow, 129626 | Russian Federation  
Tel. +7 495 663 22 11  
thommenmedical@geosoft.ru

### SINGAPORE

FONDACO Pte Ltd  
7 Kaki Bukit Road 1, #03-06  
Eunos Techno Link  
Singapore 415937 | Singapore  
Tel. +65 6392 2806  
Fax +65 6392 1296  
fondaco@fondacosg.com

### SOUTH KOREA

KMbio  
02 Ho, 129, Dongseo-daero  
Seobuk-gu, Cheonan-si  
Chungcheonnam-do  
Republic of Korea  
Tel. +82 070 3141 2875  
kmbio149@naver.com

### SPAIN/PORTUGAL

Thommen Medical Ibérica  
C/Los quintos n 1  
03350 Cox (Alicante) | Spain  
Tel. +34 96 536 10 20  
Mobile +34 606 99 78 34  
info@thommeniberica.com

### SWITZERLAND

Thommen Medical AG  
Neckarsulmstrasse 28  
2540 Grenchen | Switzerland  
Tel. +41 32 644 30 20  
Fax +41 32 644 30 25  
info@thommenmedical.ch

### TAIWAN

En-Jye International Co., Ltd.  
No. 18 | Lane 177 | Sec 3 | Chengde Rd.  
Taipei, 103 Taiwan  
Tel. +886 2 2585 1669  
Fax +886 2 2585 0892  
enjye168@gmail.com

### TURKEY

Bioport Biyolojik Maddeler A.Ş.  
Büyükdere cd. Subay evleri 9. Blok D1 Esentepe  
Şişli 34394 İstanbul | Turkey  
Tel. +90 212 2727577  
Fax +90 212 2727628  
info@bioport.com.tr  
www.bioport.com.tr

### USA/CANADA

Thommen Medical USA L.L.C.  
1375 Euclid Avenue | Suite 450  
Cleveland OH 44115 | USA  
Tel. +1 866 319 9800 (toll free)  
Fax +1 216 583 9801  
info.us@thommenmedical.com  
orders.us@thommenmedical.com