hyADENT BG

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Stand der Information: 05.2021

CE 2797

Hyaluronic acid – the biological power in dentistry

INSTRUCTIONS FOR USE

hyaDENT BG is a sterile gel based on hyaluronic acid of non-animal source.

Description

hyaDENT BG is a transparent and highly viscous gel.

hyaDENT BG is a medical device that contains highly purified hyaluronic acid of non-animal origin with a specially adjusted molecular weight.

Composition

1 ml **hyaDENT BG** contains:

Hyaluronic acid	2.0 mg
Hyaluronic acid cross-linked	16.0 mg
Sodium chloride	6.9 mg
Water for injection ad	1.0 ml

Indication

hyaDENT BG is a sterile medical device intended for single use. **hyaDENT BG** applied in the following dental disciplines is enabling wound healing:

- In treatment of gingival recession sites by coronally advanced flap (CAF) procedure.
- After tooth extraction due to a positive impact on reduction of postoperative swelling, pain and trismus.

In the environment of aesthetic dentistry **hyaDENT BG**, injected into the receded interdental papillae, leads to a reconstructive effect.

Contraindication

hyaDENT BG may not be administered to patients with

- Auto-immune disorders or under immunotherapy
- Sensitivity to hyaluronic acid
- Multiple severe allergies
- Active inflammatory or infection process
- Coagulation defects or anticoagulant therapy
- Patient with hormonal disorders

No clinical data are available on the application of **hyaDENT BG** during pregnancy or lactation or to adolescents under 18 years of age.

Adverse events

The use of **hyaDENT BG** is associated with the risk of adverse tissue reactions in very rare cases. In common with other procedures of this type the use of **hyaDENT BG** is associated with the inherent risk of infection. Papilla reconstruction may be associated with transient pressure pain due to the injection. There are no known interactions between medications and **hyaDENT BG**.

Application

One cartridge is removed from the blister and inserted into a cartridge syringe. For the application of the product on the treatment area blunt needles (23G) are recommended. For papilla reconstruction sharp needle (27G) are recommended for injection. The product is degraded completely over time.

Nature and contents of container

hyaDENT BG is supplied in a blister for single use in a 1.2 ml cartridge. The injectable content of **hyaDENT BG** has been heat sterilized. Instructions for use and labels containing the lot number and the shelf life are packaged in a carton together with the blister pack. The patient is given one of the labels to guarantee traceability of the product.

hyaDENT BG must not be re-sterilised. If the packaging is open/defective or if the blister is damaged, it must not be used as the sterility of the product is no longer guaranteed which is associated with an increased risk of infection. The package includes two cartridges á 1.2 ml.

Important information

The graduation on the syringe is intended as an orientation aid for users with regard to final volume. It does not have any measuring function, it solely provides the amount used in relation to the nominal volume of 1.2 ml.

Storage

Do not freeze or expose to extreme heat.
Store at room temperature (2 °C – 25 °C).

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BioScience GmbH