


Explantation drills for PF 3.5 – 6.0

Instructions for use THM61152


1. At a glance

This manual relates to all explantation drills for Thommen NEVO, ELEMENT, CONTACT, ONETIME and DIRECT implant variants, including the associated accessory parts. The explantation drills for NEVO, ELEMENT and CONTACT implants are listed in the product catalogue (<https://www.ifu-tm.com/THM31111>). The explantation drills for ONETIME and DIRECT are not in the product catalogue, but these components are still available. To purchase, please contact your country representative. You can find all the article numbers for the parts described here in these instructions for use.

Explantation drill, short



Explantation drill, long



Guiding cylinder for NEVO, ELEMENT and CONTACT implants, stainless steel



Components	Material
Guiding cylinder	Stainless steel
Explantation drill, short, 31 mm*	Stainless steel
Explantation drill, long, 37 mm	Stainless steel
Abutment screw	Titanium alloy
Occlusal screw for bar	Titanium alloy



*** Note:** Only use explantation drills 31.0 mm in length for the following implant lengths:

NEVO	6.5-12.5 mm
ELEMENT	6.5-12.5 mm
CONTACT	6.5-12.5 mm
CONTACT PF 6.0	9.5-11.0 mm



Important notes: The explantation drills and guiding cylinder are intended for single use only. The explantation drills and guide cylinders are supplied in a non-sterile condition and must be correspondingly processed prior to first use. See Instructions for use «THM61131 Processing». In order to ensure perfect cutting by the explantation drill, it should be used intermittently and the cutting must be continually freed from the removed bone. In order to avoid overheating the bone, the explantation drill must be used in combination with cooled physiological saline solution.

2. Notes on planning an explantation

If the explantation of an implant is not possible using the adapter for hand piece, one-piece (PF 3.5 art. no. 3.03.248, PF 4.0 art. no. 3.03.251, PF 4.5-6.0 art. no. 3.03.241), the explantation drill must be used. Specific explantation drills and guiding cylinders are required for the explantation, which are compatible with the different Thommen dental implants. The superstructure must be removed first when using these explantation drills. In order to prevent a preparation of the bone being too deep and a possible lesion of the existing nerve pathways and blood vessels, it is important to use an explantation drill which corresponds with the implant length. If the implant length is not known, the exact length of the implant and the condition of the bone can be determined using an orthopantomogram.

In preoperative planning one should consider that explantation of a CONTACT implant can be difficult because of its conical-cylindrical design, if the apical end of the implant is implanted close to the roots of adjacent teeth.



Note: If the preparation is not deep enough, because an explantation drill which is too short has been used, this can make the explantation difficult. In order to protect the patient from metal particles, which occur while the implant is being removed, a rubber dam should be used during the explantation.

Intended use

Thommen Medical instruments are used in combination with the Thommen Medical Dental Implant System in the upper and lower jaw for implant-borne tooth replacement.

Indication

Thommen Medical instruments are used in combination with the Thommen Medical Dental Implant System in partially edentulous/edentulous upper and lower jaw for restoration of masticatory function.

3. Product overview and step-by-step procedure

3.1 Explantation of NEVO and ELEMENT implants

NEVO/ELEMENT PF3,5	Rotation speed (rpm) max. 400				Work step 4a	Work step 4b
	Work step 1	Work step 2	Work step 3a	Work step 3b		
Instruments	3.03.308 Guide cylinder/ 4.03.500 Abutment screw		3.03.300 Explanation drill, long	3.03.303 Explanation drill, short	3.03.300 Explanation drill, long	3.03.303 Explanation drill, short
L 8 + 9.5 mm			This step is omitted	This step is omitted	 or 	
L 11 + 12.5 mm			This step is omitted			
L 14 mm			This step is omitted	This step is omitted		Not possible
L 17 mm (only for ELEMENT)				This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 11+12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation

NEVO/ELEMENT PF4,0	Rotation speed (rpm) max. 300				Work step 4a	Work step 4b
	Work step 1	Work step 2	Work step 3	Work step 3b		
Instruments	3.03.307 Guide cylinder/ 4.03.501 Abutment screw		3.03.304 Explanation drill, short		3.03.301 Explanation drill, long	3.03.304 Explanation drill, short
L 6.5 + 8 + 9.5 mm			This step is omitted			
L 11 + 12.5 mm						
L 14 mm			This step is omitted			Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 11 + 12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 11 + 12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation

NEVO/ELEMENT PF4.5	Rotation speed [rpm] max. 300					Work step 3	Work step 4a				Work step 4b	Work step 5a		Work step 5b
	Work step 1	Work step 2	Work step 3		Work step 4a		Work step 4a	Work step 4a	Work step 4a	Work step 4a		Work step 5a	Work step 5a	
Instruments	3.03.309 Guide cylinder/ 4.03.501 Abutment screw				3.03.301 Explantation drill, long					3.03.304 Explantation drill, short		3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	
L 6.5 + 8 + 9.5 mm					This step is omitted		This step is omitted			This step is omitted		or		
L 11 + 12.5 mm					This step is omitted		This step is omitted			*)		or		
L 14 mm					This step is omitted		This step is omitted							Not possible
L 17 mm (only for ELEMENT)					This step is omitted		*)							Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	Circular preparation of pro- jecting part at least to the di- ameter of the guide cylinder [by drilling, grinding].		*) For the 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.					*) For the 11+12.5 mm implant if the guide cylin- der is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.		Explanation	Explanation	


NEVO/ELEMENT PF5.0	Rotation speed [rpm] max. 200					Work step 3	Work step 4a		Work step 4b
	Work step 1	Work step 2	Work step 3		Work step 4a		Work step 4a	Work step 4a	
Instruments	3.03.310 Guide cylinder/ 4.03.501 Abutment screw				3.03.305 Explantation drill, short				3.03.305 Explantation drill, short
L 8.5 + 8 + 9.5 mm					This step is omitted		or		
L 11 + 12.5 mm					*)		or		
L 14 mm					This step is omitted				Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 11 + 12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.				Explanation	Explanation	


		Rotation speed (rpm) max. 200		
ELEMENT PF 6.0	Work step 1	Work step 2	Work step 3	Work step 4
Instruments	3.03.313 Guide cylinder/ 4.03.501 Abutment screw		3.03.306 Explantation drill, short	3.03.306 Explantation drill, short
L 6.5 + 8 + 9.5 mm			This step is omitted	 
L 11 + 12.5 mm			 *)	
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 11 + 12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explantation

3.2 Explantation of CONTACT implants

		Rotation speed (rpm) max. 400			
	Work step 1	Work step 2	Work step 3	Work step 4a	Work step 4b
CONTACT PF 3.5					
Instruments	3.03.308 Guide cylinder/ 4.03.500 Abutment screw		3.03.303 Explantation drill, short	3.03.300 Explantation drill, long	3.03.303 Explantation drill, short
L 11 + 12.5 mm					
L 14 mm			This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 11 + 12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation
		Rotation speed (rpm) max. 300			
	Work step 1	Work step 2	Work step 3a	Work step 3b	Work step 4a
CONTACT PF 4.0					
Instruments	3.03.307 Guide cylinder/ 4.03.501 Abutment screw		3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	3.03.301 Explantation drill, long
L 9.5 mm			This step is omitted	This step is omitted	
L 11 + 12.5 mm			This step is omitted		
L 14 mm			This step is omitted	This step is omitted	
L 17 mm				This step is omitted	
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 11+12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation

		Rotation speed (rpm) max. 300						
CONTACT PF 4,5	Work step 1	Work step 2	Work step 3	Work step 4a	Work step 4b	Work step 5a	Work step 5b	
Instruments	3.03.309 Guide cylinder/ 4.03.501 Abutment screw			3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	
L 9.5 mm				This step is omitted	This step is omitted	 or  or 		
L 11 + 12.5 mm				This step is omitted	 *)			
L 14 mm				This step is omitted	This step is omitted		Not possible	
L 17 mm				 *)	This step is omitted		Not possible	
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	Circular preparation of projecting part at least to the diameter of the guide cylinder (by drilling, grinding).	*) For the 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 11+12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation	

			Rotation speed (rpm) max. 200			
CONTACT PF 5.0	Work step 1	Work step 2	Work step 3a	Work step 3b	Work step 4a	Work step 4b
Instruments	3.03.310 Guide cylinder/ 4.03.501 Abutment screw		3.03.302 Explantation drill, long	3.03.305 Explantation drill, short	3.03.302 Explantation drill, long	3.03.305 Explantation drill, short
L 9.5 mm			This step is omitted	This step is omitted	 or 	
L 11 + 12.5 mm			This step is omitted	 *)	 or 	
L 14 mm			This step is omitted	This step is omitted		Not possible
L 17 mm			 *)	This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 11+12.5 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation

			Rotation speed (rpm) max. 200			
CONTACT PF 6.0	Work step 1	Work step 2	Work step 3	Work step 4	Work step 5a	Work step 5b
Instruments	3.03.310 Guide cylinder/ 4.03.501 Abutment screw			3.03.305 Explantation drill, short	3.03.302 Explantation drill, long	3.03.305 Explantation drill, short
L 9.5 mm + 11 mm			Important*		 or 	
L 12.5 + 14 mm			Important*	This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	Circular preparation of projecting part at least to the diameter of the guide cylinder (by drilling, grinding). * Important: CONTACT Implants with Ø shoulder 6.0 mm must be prepared up to 4 mm below the edge of the etched surface to the diameter of the guide cylinder (Ø 5.2 mm). We recommend uncovering the implant in this region in advance with a round burr.	*) For the 9.5 + 11 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation

3.3 Explantation of ONETIME implants





		Rotation speed (rpm) max. 300					
	Work step 1	Work step 2	Work step 3	Work step 4a	Work step 4b	Work step 5a	Work step 5b
ONETIME PF 5.0	3.03.309 Guide cylinder/ 4.03.501 Abutment screw			3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	3.03.301 Explantation drill, long	3.03.304 Explantation drill, short
L 8 mm				This step is omitted	This step is omitted	 or 	
L 9.5 + 11 mm				This step is omitted	 *)	 or 	
L 12.5 + 14 mm				This step is omitted	This step is omitted		Not possible
L 17 mm				 *)	This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	Circular preparation of projecting part at least to the diameter of the guide cylinder (by drilling, grinding).	*) For the 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 9.5 + 11 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	Explanation	Explanation

		Rotation speed (rpm) max. 200				
ONETIME PF 6.0	Work step 1	Work step 2	Work step 3	Work step 4	Work step 5a	Work step 5b
Instruments	3.03.310 Guide cylinder/ 4.03.501 Abutment screw			3.03.305 Explantation drill, short	3.03.302 Explantation drill, long	3.03.305 Explantation drill, short
L 8.0 mm				This step is omitted	 or 	
L 9.5 + 11 mm				 *)	 or 	
L 12.5 + 14 mm				This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	Circular preparation of projecting part at least to the diameter of the guide cylinder (by drilling, grinding).	*) For the 9,5 + 11 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.		



Notes: Guiding cylinders for the explantation of ONETIME implants are not listed in the product catalogue. But these parts are still available. To purchase, please contact your country representative.

3.4 Explantation of DIRECT implants

Rotation speed (rpm) max. 400												
			Work step 3	Work step 2	Work step 1	Work step 4a	Work step 4b	Work step 5a	Work step 5b	Work step 6a	Work step 6b	
DIRECT PF 3.5					3.03.311 Guide cylinder/4.03.506 Occlusal screw for bar	3.03.300 Explantation drill, long	3.03.303 Explantation drill, short	3.03.300 Explantation drill, long	3.03.303 Explantation drill, short	3.03.300 Explantation drill, long	3.03.303 Explantation drill, short	
L 8 mm						This step is omitted	 *)	This step is omitted	 *)	 or 	 or 	
L 11 mm						This step is omitted	 *)	This step is omitted	 *)	 or 		Not possible
L 14 mm						 *)	This step is omitted	 *)	This step is omitted			Not possible
L 17 mm						 *)	This step is omitted	 *)	This step is omitted			Not possible
	Screw on the guide cylinder using the occlusal screw for bar.	Circular preparation of projecting part at least to the diameter of the guide cylinder (by drilling, grinding).	Fold back the gingiva.			*) For the 14 + 17 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) For the 8 + 11 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the implant with the drill.	*) Shorten the implant by 1 mm for an implant of 14 mm and shorten the implant by 4 mm for an implant of 17 mm in length.		Explanation	Explanation	

		Rotation speed (rpm) max. 400						
DIRECT PF 4.2	Work step 1	Work step 2	Work step 3a	Work step 3b	Work step 4a	Work step 4b	Work step 5a	Work step 5b
Instruments	3.03.312 Guide cylinder/ 4.03.506 Occlusal screw for bar		3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	3.03.301 Explantation drill, long	3.03.304 Explantation drill, short	3.03.301 Explantation drill, long	3.03.304 Explantation drill, short
L 8 mm			This step is omitted	 *)	This step is omitted	 *)	 or 	 or 
L 11 mm			This step is omitted	 *)	This step is omitted	 *)	 or 	 or 
L 14 mm			 *)	This step is omitted	 *)	This step is omitted		Not possible
L 17 mm			 *)	This step is omitted	 *)	This step is omitted		Not possible
	Screw on the guide cylinder using the abutment screw.	Fold back the gingiva.	*) For the 14 + 17 mm im- plant if the guide cylinder is bumped against, this guide cylinder is removed in order to be able to reach the tip of the im- plant with the drill.	*) For the 8 + 11 mm implant if the guide cylinder is bumped against, this guide cylinder is removed in or- der to be able to reach the tip of the implant with the drill.	*) Shorten the implant by 1 mm for an implant of 14 mm and shorten the im- plant by 4 mm for an im- plant of 17 mm in length.	*) Shorten the implant by 1 mm for an implant of 8 mm and shorten the implant by 4 mm for an implant of 11 mm in length.	Explanation	Explanation



Notes: Guiding cylinders for the explantation of DIRECT implants are not listed in the product catalogue. But these parts are still available. To purchase, please contact your country representative.

General notes

THOMMEN IMPLANT SYSTEM

	Manufacturer: Thommen Medical AG Neckarsulmstrasse 28 2540 Grenchen, Switzerland www.thommenmedical.com
	Batch code
	Use by date
	Date of manufacture
	Sterilized using irradiation
	Simple sterile barrier system
	Simple sterile barrier system with protective packaging on the inside
	Not more than 20 reprocessing cycles
	European representative
	Temperature limitation
	Do not re-use
	Non-sterile
	Caution
	Article number
	Conformity symbol in accordance with Directive No: MDD 93/42/EEC or Regulation (EU) No: 2017/745 MDR (see corresponding conformity declaration)
	Consult instructions for use
	Do not resterilize
	Do not use if package is damaged
	Atmospheric pressure limitation
	Manufacturer
	Country of manufacture
	Importer
	Dealers
	Keep away from sunlight
	May only be sold to and prescribed by physicians (USA)
	Medical device
	Single product code

COLORED WARNING STICKER

Application was changed – follow the directions in the corresponding instructions for use.

NEW HANDLING

New design – the application has not been changed.

NEW DESIGN

PRODUCT INFORMATION The information in this document describes the application of the Thommen Medical implant system. This information is available in electronic form online at: www.ifu-tm.com. The responsible country representative or distributor for Thommen Medical AG is available to provide technical advice.

COLOR CODE Each implant platform diameter has a color code, which can be found on all implant packagings, on the impression items and on most diameter-specific instruments.

TRACEABILITY

In order to ensure the traceability of the implantable products as well as the manufacturer, product type and product dimensions also for a later prosthetic re-restoration, each product package comes with three patient labels. These labels should be used in the practice for documentation and for the implant passport.

Brown	=	PF 3.0
Yellow	=	PF 3.5
Green	=	PF 4.0
Blue	=	PF 4.5
Grey	=	PF 5.0
Purple	=	PF 6.0

AVAILABILITY Not all of the Thommen Medical products mentioned in these instructions for use are available in all countries. The responsible country representative or distributor for Thommen Medical AG can provide information about availability of Thommen Medical products for the country in question.

GENERAL RESTRICTIONS OF USE Restorations with cantilevers to individual implants are not recommended. Individual restorations with angled abutments should not be used in regions with high mechanical stress. For implants with a small diameter (PF 3.0 and 3.5), the prosthetic restoration should be constructed in such a way that a large bending moment does not occur.

CONTRAINDICATION The Thommen Medical products may not be used on patients who are known to have allergies to the corresponding materials.

POSSIBLE COMPLICATIONS A stressed loading of the implant or abutment over and above its functional capacity can lead to excessive bone loss or fracture of the implant or restoration. The clinician must supervise the occlusion and functional loading of the prosthetic supraconstruction very carefully.

SIDE EFFECTS The patient should be informed about the possible side effects, interactions, precautionary measures and complications associated with Thommen Medical products. Potential complications can occur immediately after insertion of dental implants: Temporary symptoms: swelling, difficulties with speaking, gum inflammations, pain. Longer lasting symptoms: chronic pain connected with the dental implant, localized or systemic infections, dysesthesia, loss of alveolar ridge (upper and lower jaw), oronasal or oronasal fistulas, irreversible damage to neighboring teeth, esthetic problems, nerve damage, hyperplasia. Other possible side effect: implant loss.

WARNINGS All Thommen Medical products that come into effect inside the oral cavity must be protected against aspiration. Thommen Medical products have not been tested for safety and compatibility in an MR environment. Thommen Medical products have not been tested for heating or migration in the MR environment. The safety of Thommen Medical products in the MR environment is unknown. Magnetic resonance tomographic examinations of patients, who have been treated with Thommen Medical products, may result in patient injuries.

RESPONSIBILITY/LIABILITY As a part of an overall scheme, Thommen Medical products may be used only with the related original components and instruments in accordance with the instructions for use provided by Thommen Medical. The use of non-system parts may compromise the performance of Thommen Medical products and lead to failures. Users must have appropriate knowledge and information about the handling of Thommen Medical products in order to use the products safely and correctly. The user is obliged to use the Thommen Medical products according to the instructions for use and to check whether the product is suitable for the individual patient situation. The use of Thommen Medical products is the responsibility of the user, as such, beyond the control of Thommen Medical AG. We refuse to accept any responsibility or liability for any damage due to incorrect utilization of the product. Products labeled "Do not re-use" may not be refurbished and/or reused. The refurbishment and/

or reuse of these products can affect their function (e.g. fitting and/ or cutting properties) as well as their safe use (e.g. risk of infection, disease transmission, fading of the laser or color marks, corrosion). Detailed information about the possible consequences, which may result from incorrect use, is available from the responsible country representative or distributor of Thommen Medical AG. All serious incidents which have occurred in connection with the product must be reported to the manufacturer and to the local competent authority in the country where the user is resident.

GUARANTEE The comprehensive guarantees can be found in the country-specific guarantee leaflets.

TRANSPORT AND STORAGE Please note the specifications on the labels and instructions for use regarding transportation, storage and handling. If the packaging is damaged, the products must not be used; a visual inspection is necessary. Under no circumstances may Thommen Medical products be used beyond the expiry date, as proper functioning or sterility of sterile packaged products cannot be guaranteed by the manufacturer.

APPLICATION The following descriptions are not intended as comprehensive for the immediate use of the Thommen Medical Implant System. Training by a specialist experienced in the use of this system is recommended.

GUARANTEE OF STERILITY In general, products of the Thommen Implant System supplied in sterile packaging must not be re-sterilized. Sterile-packed products, whose packaging is damaged, must not be used under any circumstances. Sterile-supplied products, which have not been used for the surgical operation, whose packaging has been opened are considered as having been used and must not be used thereafter. In the event of re-sterilization, proper function and the sterility cannot be guaranteed by the manufacturer. The products intended for single use must never be reprocessed, sterilized or reused and must be disposed of safely and properly after use in compliance with all applicable legal and regulatory requirements. Reusable products must be reprocessed according to the instructions for use and, if used on patients, sterilized. They must be checked for their integrity before each use. Any damage (for example, scratches, cracks, nicks, dents), as well as bent parts, means that they must not be used any longer. The number of reprocessing cycles is limited and must be monitored. If the number of cycles is exceeded, proper function and sterility of the product are not guaranteed by the manufacturer anymore.

DISPOSAL In the case of cutting products, there is always a risk of injury, therefore the products must be disposed of safely and properly after use, observing all applicable legal and regulatory requirements. For products and their accessories, which have been used on a patient, there is a risk of an infection. Our products are designed and produced so that they can be disposed of safely and correctly after use in compliance with all valid legal and regulatory requirements.

COPYRIGHT®/REGISTERED® SPI®, INICELL® and APLIQUIQ® are registered trademarks of Thommen Medical AG. Publication or reproduction is permitted only with the written consent of Thommen Medical AG. LOCATOR® is a registered trademark of Zest Anchors Inc., CA, USA. Dynamic Tibase® is a registered trademark of Talladium España S.L., Lleida, Spain.

VALIDITY Thommen Medical AG. All previous versions lose their validity with the publication of this instruction for use.



HEADQUARTERS

Thommen Medical AG
Neckarsulmstrasse 28
2540 Grenchen | Switzerland
Tel. +41 61 965 90 20
Fax +41 61 965 90 21
info@thommenmedical.com



SUBSIDIARIES/NATIONAL DISTRIBUTORS

AUSTRIA

Thommen Medical Austria GmbH
Mühlgasse 3
2322 Zwölfaxing | Austria
Tel. +43 660 2011953
info@thommenmedical.at

BENELUX

Thommen Medical Benelux B.V.
Dierenriem 1
3738 TP Maartensdijk | Netherlands
Tel. +31 30 68 68 468
Info.benelux@thommenmedical.nl

CHINA

Shanghai Yujing Trading Co., Ltd.
Room 1112 | Block A | SOHO Zhongshan Plaza
No. 1055 West Zhongshan Rd | Changning District
200051 Shanghai | P.R. China
Tel. +86 21 80121250
Fax +86 21 62175264
sy@shanghaiyujing.com

CZECH REPUBLIC

C. Witt Dental spol. s r.o.
Cihlářská 643/19
602 00 Brno
Tel. +420 739 043 449
helen.novak@cwittdental.cz

FINLAND

Vector Laboratories Oy
Engelinaukio 8 B
00150 Helsinki | Finland
Tel. +358 400 940 700
labs@vektor.fi

FRANCE

Thommen Medical France
Le PARK, Bâtiment B, 1 Rue Charles Cordier
77164 Ferrières-en-Brie | France
Tel. +33 1 83 64 06 35
infos@thommenmedical.fr
commande@thommenmedical.fr



GERMANY

Thommen Medical Deutschland GmbH
Am Rathaus 2
79576 Weil am Rhein | Germany
Tel. +49 7621 422 58 30
Fax +49 7621 422 58 41
info@thommenmedical.de

HONG KONG

Shengyuan (Hong Kong) Int. Trade Co. Ltd.
Level 13, 68 Yee Wo Street
Causeway Bay | Hong Kong
Tel. +852 530 876 41

HUNGARY

Mori-Dent Kft.
Huszt utca 9 | 1147 Budapest
Hungary
Tel. +361 220 5443
morident@morident.hu
www.morident.hu

ITALY

Dental Trey S.r.l.
Via Partisani, 3
47016 Fiumana | Predappio (FC) | Italy
Tel. +39 0543 929111
Fax +39 0543 940659
implantologia@dental Trey.it
www.dental Trey.it

JAPAN

J. Morita Corporation
3-33-18, Tarumi-cho
Suita | Osaka 564-8650 | Japan
Tel. +81 6 6384 6921
Fax +81 6 6384 6746
www.morita.com

JORDAN

Alzynoor for Medical Supplies
Complex No. 243
Barakat Center office 406, JO-Amman 11953 Jordan
Tel. +962 655 39020
alzynoor.medical.jo@gmail.com

KUWAIT

Logident
Yasmin Complex, 1st floor, Office B1
Ibn Khaldoun, St Hawally | Kuwait
Tel. +965 5113 8438
omar.alameddin@logidentco.com

LITHUANIA/LATVIA

LT projects, UAB
Europos pr. 121 | Kaunas
LT-46339 Lithuania
Tel. +370 65 771550
info@fordentist.lt
www.fordentist.lt

NORWAY

Novus Dental AS
Måkeveien 6
1679 Kråkerøy. | Norway
post@novusdental.no
www.novusdental.no

OMAN

Al Inshirah Street, Madinat Sultan Qaboos
Building No.185, Way No.1947
P.O. Box 1181, P.C. 112.
Ruwi, Muscat,
Sultanate Of Oman
Tel. +968 24602771
sam@almadinagroup.net

POLAND

C.WITT DENTAL Sp. z o. o.
Ul. Granitowa 10
87-100 Toruń | NIP 951-15-08-371 | Poland
Tel. +48 56 623 61 23
biuro@cwittdental.pl
www.cwittdental.pl

QATAR

Orient gate company
Al - Nasr, Area 39, Street 840, Building No: 53
P.O. Box 45393 | Doha Qatar
Mobile +974 50389808
lyaseen@orientgateqa.com

REPUBLIC OF CROATIA

Futura Dental d.o.o.
Ulica grada Vukovara 271 | Chromosov toranj
10 000 Zagreb | Republic of Croatia
Tel. +385 91 6814 860
info@futura-dental.hr
www.futura-dental.hr

RUSSIAN FEDERATION

Geosoft Surgery LLC.
Vasnetsova Lane, House 7, Office 206
129090 Moscow | Russian Federation
Tel. +7 495 663 22 11
thommenmedical@geosoft.ru
www.geosoft.ru

SAUDI ARABIA

INNOVATION MEDICAL COMPANY
7324 King Abdulaziz Road
Ar-Rabi District
13315 Riyadh | Kingdom of Saudi Arabia
Tel. +966 50 926 5010
contact@innovation-medical.com

SINGAPORE

FONDACO Pte Ltd
7 Kaki Bukit Road 1, #03-06
Eunos Techno Link
Singapore 415937 | Singapore
Tel. +65 6392 2806
Fax +65 6392 1296
fondaco@fondacosg.com

SOUTH KOREA

KMbio
02 Ho, 129, Dongseo-daero
Seobuk-gu, Cheonan-si
Chungcheongnam-do
Republic of Korea
Tel. +82 070 3141 2875
kmbio149@naver.com

SPAIN/PORTUGAL

Thommen Medical Ibérica
C/Los quintos n 1
03350 Cox (Alicante) | Spain
Tel. +34 96 536 10 20
Mobile +34 606 99 78 34
info@thommeniberica.com

SWEDEN

Erik Söderberg Group AB
Björnkärsvägen 112
187 42 Täby | Sweden
Tel. +46 70 567 00 67
erik@thommenmedical.se
www.thommenmedical.se

SWITZERLAND

Thommen Medical AG
Neckarsulmstrasse 28
2540 Grenchen | Switzerland
Tel. +41 32 644 30 20
Fax +41 32 644 30 25
info@thommenmedical.ch

TAIWAN

En-Jye International Co., Ltd.
No. 18, Lane 177, Sec 3, Chengde Rd.
Taipei, 103 Taiwan | Taiwan
Tel. +886 2 2585 1669
Fax +886 2 2585 0892
enjye168@gmail.com

TURKEY

Bioport Biyolojik Maddeler A.Ş.
Büyükdere cd. Subay evleri 9. Blok D1 Esentepe
Şişli 34394 İstanbul | Turkey
Tel. +90 212 2727577
Fax +90 212 2727628
info@bioport.com.tr
www.bioport.com.tr

UKRAINE

BIG TIME GROUP LTD
Kyivska 5/1A | 46016 Ternopil
Ukraine
Tel. +380 67 54 60 147
btgukraine@gmail.com
bigtimegroup.com.ua

UNITED ARAB EMIRATES

City Pharmacy Co WLL
Taweelah, KIZAD Sector KHIA 9 | Plots 68 & 69
P.O.Box 2098, Abu Dhabi | UAE
Tel. +971 507 407 047 or +971 5056 20568
faisal.m@citypharmacy.com
moawia.zaharna@citypharmacy.com

USA/CANADA

Thommen Medical USA L.L.C.
1375 Euclid Avenue | Suite 450
Cleveland OH 44115 | USA
Tel. +1 866 319 9800 (toll free)
Fax +1 216 583 9801
info@thommenmedical.us
orders@thommenmedical.us