


pluser

ERBIUM LASER TIP CATALOGUE

Each pack of Tip contains 5 pieces

		SOFT TISSUE	TOOTH	BONE	AESTHETIC	PERIO	ENDO
 Specific tip for bone surgery, periodontics and peri-implantitis.	CODE LAEH4008.5/A DIAMETER 400μ LENGTH 8 mm			✓ BONE		✓ PERIO	
 Specific tip for soft tissue and bone surgery. Cased.	CODE LAEH4008.5/A1 DIAMETER 400μ LENGTH 8 mm	✓ SOFT TISSUE		✓ BONE			
 Specific tip for bone surgery, periodontics and peri-implantitis.	CODE LAEH4012.5/A DIAMETER 400μ LENGTH 12 mm			✓ BONE		✓ PERIO	
 Specific tip for soft tissue and bone surgery. Cased.	CODE LAEH4012.5/A1 DIAMETER 400μ LENGTH 12 mm	✓ SOFT TISSUE		✓ BONE			

Each pack of Tip contains 5 pieces		SOFT TISSUE	TOOTH	BONE	AESTHETIC	PERIO	ENDO
 Specific tip for bone surgery and posterior sectors.	CODE LAEH6004.5/A DIAMETER 600μ LENGTH 4 mm			✓ BONE			
 Specific tip for soft tissue and bone surgery, periodontics and peri-implantitis.	CODE LAEH6008.5/A DIAMETER 600μ LENGTH 8 mm	✓ SOFT TISSUE		✓ BONE		✓ PERIO	
 Specific tip for soft tissue and bone surgery. Cased.	CODE LAEH6008.5/A1 DIAMETER 600μ LENGTH 8 mm	✓ SOFT TISSUE		✓ BONE			
 Specific tip for bone surgery, periodontics and peri-implantitis in anterior sectors.	CODE LAEH6012.5/A DIAMETER 600μ LENGTH 12 mm			✓ BONE		✓ PERIO	
 Specific tip for soft tissue and bone surgery. Cased.	CODE LAEH6012.5/A1 DIAMETER 600μ LENGTH 12 mm	✓ SOFT TISSUE		✓ BONE			
 Specific tip for bone surgery and periodontics. Cased.	CODE LAEH6020.5/A1 DIAMETER 600μ LENGTH 20 mm			✓ BONE		✓ PERIO	
 Specific tip for periodontics and peri-implantitis. Cased	CODE LAEH6025.5/A1 DIAMETER 600μ LENGTH 25 mm					✓ PERIO	

Each pack of Tip contains 5 pieces		SOFT TISSUE	TOOTH	BONE	AESTHETIC	PERIO	ENDO
	CODE LAEH8004.5/A DIAMETER 800μ LENGTH 4 mm Specific tip for conservative, bone surgery and aesthetic.		✓ TOOTH	✓ BONE	✓ AESTHETIC		
	CODE LAEH8008.5/A DIAMETER 800μ LENGTH 8 mm Tip for endodontics through BLED technique (Bubble Laser Endodontics Decontamination).						✓ ENDO
	CODE LAEH8008.5/A1 DIAMETER 800μ LENGTH 8 mm Specific tip for soft and bone surgery, conservative and aesthetic. Cased.	✓ SOFT TISSUE	✓ TOOTH	✓ BONE	✓ AESTHETIC		
	CODE LAEH8012.5/A DIAMETER 800μ LENGTH 12 mm Specific tip for bone surgery and conservative.		✓ TOOTH	✓ BONE			
	CODE LAEH8012.5/A1 DIAMETER 800μ LENGTH 12 mm Specific tip for bone surgery and conservative. Cased.		✓ TOOTH	✓ BONE			
	CODE LAEH8020.5/A1 DIAMETER 800μ LENGTH 20 mm Specific tip for periodontics and peri-implantitis. Cased.					✓ PERIO	



Each pack of Tip contains 5 pieces

		SOFT TISSUE	TOOTH	BONE	AESTHETIC	PERIO	ENDO
 Specific tip for bone surgery. Chisel tip for thin cutting.	CODE LAEH8081.5/A DIAMETER 800x300μ LENGTH 8 mm			 BONE			
 Specific tip for bone surgery. Chisel tip for thin cutting. Cased.	CODE LAEH8081.5/A1 DIAMETER 800x300μ LENGTH 8 mm			 BONE			

LASER THERAPY

ALL ROUND

Designed to be precise on even the most complex tissues, the **tips** provided with Pluser help to achieve fantastic results, working not in contact and allowing patients to avoid all unpleasant vibrations from traditional rotating instruments.

The tips, which are identified by colour code, have specific lengths and diameters to ensure **the most efficient laser emission** according to the treatment to be carried out, with all the benefits of the best and minimally invasive treatment for your patients.

Coloured tips according to treatment

Sapphire tips for top performance

Quick and easy insertion into the handpiece



LAMBDA SpA
Via dell'Impresa, 1 | 36040 Brendola (VI) Italy
T +39 0444.349165 | F +39 0444.349954
info@lambdaspa.com | lambdaspa.com



AMDPI196.0



INFO

• doctor-smile.com • info@doctor-smile.com

All LAMBDA S.p.A. laser equipment
has EC certification pursuant to EEC Directives
for medical equipment: ECC Directive 93/42

